KOMO

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Center for Telemedicin, Region Midtjylland
Olof Palmes Allé 36 8200 Århus

1 Intended use

Healthcare professionals use the system for setup of questionnaires to answer and measurements to be performed by patients and for reviewing questionnaire responses and measurement results reported by patients. Healthcare professionals attach threshold values to questionnaires and measurement regimes and use the threshold values for decision support when viewing questionnaire responses and measurement results. These threshold values are used to triage measurement results and questionnaire responses into three categories: red (pathological answer or measurement), yellow (deviant answer or measurement), and green (normal answer or measurement).

When presenting a list of patient responses awaiting review and feedback, this categorization, if present, is indicated alongside the responses and the categorisation is furthermore used to sort the list of patients and their responses. The intended use requires that **all** patients **eventually** get their review and feedback **in a timely manner**, and that **acute matters are kept separate** from the use of KOMO ver. 2.0.0. This implies that even a "red priority" measurement or questionnaire response will not be guaranteed a fast response, only that it may be prioritized higher than the yellow and green categories. The patients are informed off these terms.

1.1 Medical indication

The indication is to monitor the well-being of the patient during treatment with intravenous antibiotic where the condition of these patients is non-time-critical. It is up to medical assessment by the treating physician whether the patient is time or non-time-critical.

1.2 Patient population

The intended patient population is treatment for patients in Central Denmark Region who need intravenous antibiotics, which can be handled by the patient themselves or a caregiver in the patient's home.

The patients are trained in how to use KOMO ver. 2.0.0 and to administrate the medication by the healthcare professional before the patient are inrolled in the treatment. During the training the patients also receive instructions that they need to contact the healthcare professionals by phone if their condition gets time-critical.

1.3 User profile

KOMO ver. 2.0.0 is intended to be used by healthcare professionals and patients for monitoring the well-being of the patient, during treatment with intravenous antibiotic handled by the patient themselves or a caregiver.

1.4 Environment of use (patients)

The IV Citizens Client of KOMO ver. 2.0.0 is a module that supports home self monitoring of patients. Patients use the system for performing measurements and/or responding to questionnaires at home, typically as an alternative to outpatient clinic visits.

The IV Citizens Client is typically used on a mobile device through which questionnaires can be answered and measurement results can be entered. Furthermore, patients can view feedback on whether a healthcare professional has reviewed the given responses and results or not and patients can also view graphs of own historical measurement results. Measurements are taken using an external measurement device. Results from such devices can be manually entered on the mobile device or the measurement device can be connected to the mobile device, transferring measurement results through wired or wireless technology (e.g. Bluetooth).

1.5 Environment of use (healthcare professionals)

The IV Employee Client of KOMO ver. 2.0.0 is intended for use in a system context with the IV Citizens Client. In the IV Employee Client the healthcare professionals use questionnaire responses and measurement result from the IV Citizens Client for monitoring the well-being of the patient, during treatment with intravenous antibiotic handled by the patient himself or a caregiver at the patients home.

The IV Employee Client is used on a computer or laptop through which questionnaire responses and measurements can be viewed and. Furthermore healthcare professionals can view graphs of the patients historical measurement results.

1.6 Contraindications

The module does not provide any real time alerts and is not intended to provide automated treatment decisions or diagnosis. The triaging thresholds shall not be used for filtering/hiding low priority ("normal") observations and questionnaire responses.

Some additional reservations:

- Not for time-critical real-time monitoring of patients' acute status
- Not for use in any time critical or emergency situations
- Not for alarms or direct patient monitoring
- Not for automated treatment or diagnosis decisions

2 Conclusion

Based on the above description, it is hereby declared that KOMO ver. 2.0.0 qualifies as a medical device, cf. EU Regulation 2017/745 on medical devices (MDR).

3 References

- REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017

midt
regionmidtjylland